



PRODUCT CATALOG

2024 STERI-OSS CATALOG -Rev03






CONTENTS

FIXTURE DESIGN
BONE LEVEL SYSTEM
TISSUE LEVEL SYSTEM
SURGICAL KIT & INSTRUMENT
GBR PRODUCT
SURGICAL INSTRUMENT

Fixture Design.



Connection


- Submerged type implant with an internal hex 11° taper connection structure
- Optimum screw thread design for SLA surface to achieve consistent roughness
- Taper body design with superior initial stability
- Highest initial stability in soft bone by using upper-section thread and body design
- Powerful self-threading effect using cutting edge design
- Acquires the initial stability needed in immediate loading even in soft bone

Cutting Edge Design

The cutting edge which applied to fixture designed to have optimum torque at the all the bone conditions. It is applied and divided by length and size. Added cutting edge to the body thread design could make convenient work for the modification of path and self tapping.

S.L.A Surface

Surface of fixture are excellent in morphology and roughness, increased 50% of the surface area than RBM. The SLA surface and the entire taper shape of Steri-oss fixture could get enough initial fix.



PACKAGE DESIGN



01

Open the package box by pushing the push part in the front side



02

Separate the cover by turning the orange color Outer Cap



03

Take out the Ampule



04

It is more conveniently could open if you rotate the transparent cover with the product upward



05

Connect the Handpiece connector with Fixture



06

Should check the upper part of fixture meet with low part of Marking Line



07


Check the Cover Screw by open the white silicone cap in the opposite direction





08

Connect the Cover screw with Fixture using 1.2 Hex driver.

STERIOSS IMPLANT BL(Bone Level) SYSTEM Reference Code.



- Surface treatment : SLA
- Internal Hex : 2.5mm
- Taper connector : 11°
- Female screw : M2.0



| Length Ø | F3.5(Ø3.8) | F4.0 (Ø4.2) | F4.5(Ø4.6) | Length Ø | F5.0(Ø5.1) | F6.0 (Ø6.0) |
|----------|------------|-------------|------------|----------|------------|-------------|
| 5 | | | SD4505S | 4 | SD5004S | SD6004S |
| 6 | | | SD4506S | 5 | SD5005S | SD6005S |
| 7 | SD3507 | SD4007 | SD4507 | 6 | SD5006 | SD6006 |
| 8 | SD3508 | SD4008 | SD4508 | 7 | SD5007 | SD6007 |
| 9 | SD3509 | SD4009 | SD4509 | 8 | SD5008 | SD6008 |
| 10 | SD3510 | SD4010 | SD4510 | 9 | SD5009 | SD6009 |
| 11 | SD3511 | SD4011 | SD4511 | 10 | SD5010 | SD6010 |
| 12 | SD3512 | SD4012 | SD4512 | 11 | SD5011 | SD6011 |
| 13 | SD3513 | SD4013 | SD4513 | 12 | SD5012 | SD6012 |
| | | | | 13 | SD5013 | SD6013 |

Cover Screw



- Long size is used in deep implant insertion
- 1.2mm Hex driver

STERIOSS IMPLANT Reference Code.


Healing Abutment



- 1.2mm Hex driver
- Recommended tightening torque : 5~8Ncm

| Height Ø | 4.0 | 4.5 | 5.0 | 5.5 | 6.0 | 6.5 |
|-----------------------|-----------|------------|------------|-----------|------------|-----------|
| 3 | AM2HAR403 | AM2HAR453 | AM2HAR503 | AM2HAR553 | AM2HAR603 | AM2HAR653 |
| 4 | AM2HAR404 | AM2HAR454 | AM2HAR504 | AM2HAR554 | AM2HAR604 | AM2HAR654 |
| 5 | AM2HAR405 | AM2HAR455 | AM2HAR505 | AM2HAR555 | AM2HAR605 | AM2HAR655 |
| 6 | AM2HAR406 | AM2HAR456 | AM2HAR506 | AM2HAR556 | AM2HAR606 | AM2HAR656 |
| 7 | AM2HAR407 | AM2HAR457 | AM2HAR507 | AM2HAR557 | AM2HAR607 | AM2HAR657 |
| 8 | AM2HAR408 | AM2HAR458 | AM2HAR508 | AM2HAR558 | AM2HAR608 | AM2HAR658 |
| 9 | AM2HAR409 | AM2HAR459 | AM2HAR509 | AM2HAR559 | AM2HAR609 | AM2HAR659 |
| 10 | | AM2HAR4510 | AM2HAR5010 | | AM2HAR6010 | |
| 11 | | AM2HAR4511 | AM2HAR5011 | | AM2HAR6011 | |


Connection Guide



All abutment except the cover screw could use to all Fixture in common.
But the over 4.0mm connection are 0.5mm higher than 3.5mm connection
For example if use 4.0mm height healing abutment to Ø4.0mm fixture the height to upper part are 4mm but it shall be 4.5mm if it us to Ø3.5mm Fixture


PROSTHETICS FLOW DIAGRAMS.

Screw & Cement Retained Restoration



Reference Code.

Solid Abutment




- Used in producing cement-retained prosthetics
- Gold coloring for improved aesthetics
- 1.2mm Hex driver
- Recommended tightening torque : 30Ncm

| Ø | | 4.0 | 4.5 | 5.0 | 6.0 |
|--------|--------|--------------|--------------|--------------|--------------|
| PH 4 | Cuff 1 | AM2SAR401040 | AM2SAR451040 | AM2SAR501040 | AM2SAR601040 |
| | 2 | AM2SAR402040 | AM2SAR452040 | AM2SAR502040 | AM2SAR602040 |
| | 3 | AM2SAR403040 | AM2SAR453040 | AM2SAR503040 | AM2SAR603040 |
| | 4 | AM2SAR404040 | AM2SAR454040 | AM2SAR504040 | AM2SAR604040 |
| | 5 | AM2SAR405040 | AM2SAR455040 | AM2SAR505040 | AM2SAR605040 |
| | 6 | AM2SAR406040 | AM2SAR456040 | AM2SAR506040 | AM2SAR606040 |
| PH 5.5 | 1 | AM2SAR401055 | AM2SAR451055 | AM2SAR501055 | AM2SAR601055 |
| | 2 | AM2SAR402055 | AM2SAR452055 | AM2SAR502055 | AM2SAR602055 |
| | 3 | AM2SAR403055 | AM2SAR453055 | AM2SAR503055 | AM2SAR603055 |
| | 4 | AM2SAR404055 | AM2SAR454055 | AM2SAR504055 | AM2SAR604055 |
| | 5 | AM2SAR405055 | AM2SAR455055 | AM2SAR505055 | AM2SAR605055 |
| | 6 | AM2SAR406055 | AM2SAR456055 | AM2SAR506055 | AM2SAR606055 |
| PH 7 | 1 | AM2SAR401070 | AM2SAR451070 | AM2SAR501070 | AM2SAR601070 |
| | 2 | AM2SAR402070 | AM2SAR452070 | AM2SAR502070 | AM2SAR602070 |
| | 3 | AM2SAR403070 | AM2SAR453070 | AM2SAR503070 | AM2SAR603070 |
| | 4 | AM2SAR404070 | AM2SAR454070 | AM2SAR504070 | AM2SAR604070 |
| | 5 | AM2SAR405070 | AM2SAR455070 | AM2SAR505070 | AM2SAR605070 |
| | 6 | AM2SAR406070 | AM2SAR456070 | AM2SAR506070 | AM2SAR606070 |

STERIOSS IMPLANT Reference Code.

Straight Abutment (HEX)




- Used in producing cement-retained prosthetics
- Gold coloring for improved aesthetics
- 1.2mm Hex driver
- Recommended tightening torque : 30Ncm

| | | Ø | 4.0 | 4.5 | 5.0 | 5.5 | 6.0 | 6.5 |
|-----------|------|---|---------------|---------------|---------------|---------------|---------------|---------------|
| PH 4 | Cuff | | | | | | | |
| | 1 | | AM2STR401040H | AM2STR451040H | AM2STR501040H | AM2STR551040H | AM2STR601040H | AM2STR651040H |
| | 2 | | AM2STR402040H | AM2STR452040H | AM2STR502040H | AM2STR552040H | AM2STR602040H | AM2STR652040H |
| | 3 | | AM2STR403040H | AM2STR453040H | AM2STR503040H | AM2STR553040H | AM2STR603040H | AM2STR653040H |
| | 4 | | AM2STR403040H | AM2STR454040H | AM2STR504040H | AM2STR554040H | AM2STR604040H | AM2STR654040H |
| | 5 | | AM2STR405040H | AM2STR455040H | AM2STR505040H | AM2STR555040H | AM2STR605040H | AM2STR655040H |
| | 6 | | AM2STR406040H | AM2STR456040H | AM2STR506040H | AM2STR556040H | AM2STR606040H | AM2STR656040H |
| PH 5.5 | 1 | | AM2STR401055H | AM2STR451055H | AM2STR501055H | AM2STR551055H | AM2STR601055H | AM2STR651055H |
| | 2 | | AM2STR402055H | AM2STR452055H | AM2STR502055H | AM2STR552055H | AM2STR602055H | AM2STR652055H |
| | 3 | | AM2STR403055H | AM2STR453055H | AM2STR503055H | AM2STR553055H | AM2STR603055H | AM2STR653055H |
| | 4 | | AM2STR403055H | AM2STR454055H | AM2STR504055H | AM2STR554055H | AM2STR604055H | AM2STR654055H |
| | 5 | | AM2STR405055H | AM2STR455055H | AM2STR505055H | AM2STR555055H | AM2STR605055H | AM2STR655055H |
| | 6 | | AM2STR406055H | AM2STR456055H | AM2STR506055H | AM2STR556055H | AM2STR606055H | AM2STR656055H |
| PH 7 | 1 | | AM2STR401070H | AM2STR451070H | AM2STR501070H | AM2STR551070H | AM2STR601070H | AM2STR651070H |
| | 2 | | AM2STR402070H | AM2STR452070H | AM2STR502070H | AM2STR552070H | AM2STR602070H | AM2STR652070H |
| | 3 | | AM2STR403070H | AM2STR453070H | AM2STR503070H | AM2STR553070H | AM2STR603070H | AM2STR653070H |
| | 4 | | AM2STR403070H | AM2STR454070H | AM2STR504070H | AM2STR554070H | AM2STR604070H | AM2STR654070H |
| | 5 | | AM2STR405070H | AM2STR455070H | AM2STR505070H | AM2STR555070H | AM2STR605070H | AM2STR655070H |
| | 6 | | AM2STR406070H | AM2STR456070H | AM2STR506070H | AM2STR556070H | AM2STR606070H | AM2STR656070H |
| | 7 | | AM2STR407070H | AM2STR457070H | AM2STR507070H | AM2STR557070H | AM2STR607070H | AM2STR657070H |

STERIOSS IMPLANT Reference Code.

Straight Abutment (NON-HEX)




- Used in producing cement-retained prosthetics
- Gold coloring for improved aesthetics
- 1.2mm Hex driver
- Recommended tightening torque : 30Ncm

| Cuff | | Ø | 4.0 | 4.5 | 5.0 | 5.5 | 6.0 | 6.5 |
|--------|---|---|---------------|---------------|---------------|---------------|---------------|---------------|
| PH 4 | 1 | | AM2STR401040N | AM2STR451040N | AM2STR501040N | AM2STR551040N | AM2STR601040N | AM2STR651040N |
| | 2 | | AM2STR402040N | AM2STR452040N | AM2STR502040N | AM2STR552040N | AM2STR602040N | AM2STR652040N |
| | 3 | | AM2STR403040N | AM2STR453040N | AM2STR503040N | AM2STR553040N | AM2STR603040N | AM2STR653040N |
| | 4 | | AM2STR403040N | AM2STR454040N | AM2STR504040N | AM2STR554040N | AM2STR604040N | AM2STR654040N |
| | 5 | | AM2STR405040N | AM2STR455040N | AM2STR505040N | AM2STR555040N | AM2STR605040N | AM2STR655040N |
| | 6 | | AM2STR406040N | AM2STR456040N | AM2STR506040N | AM2STR556040N | AM2STR606040N | AM2STR656040N |
| PH 5.5 | 1 | | AM2STR401055N | AM2STR451055N | AM2STR501055N | AM2STR551055N | AM2STR601055N | AM2STR651055N |
| | 2 | | AM2STR402055N | AM2STR452055N | AM2STR502055N | AM2STR552055N | AM2STR602055N | AM2STR652055N |
| | 3 | | AM2STR403055N | AM2STR453055N | AM2STR503055N | AM2STR553055N | AM2STR603055N | AM2STR653055N |
| | 4 | | AM2STR403055N | AM2STR454055N | AM2STR504055N | AM2STR554055N | AM2STR604055N | AM2STR654055N |
| | 5 | | AM2STR405055N | AM2STR455055N | AM2STR505055N | AM2STR555055N | AM2STR605055N | AM2STR655055N |
| | 6 | | AM2STR406055N | AM2STR456055N | AM2STR506055N | AM2STR556055N | AM2STR606055N | AM2STR656055N |
| PH 7 | 1 | | AM2STR401070N | AM2STR451070N | AM2STR501070N | AM2STR551070N | AM2STR601070N | AM2STR651070N |
| | 2 | | AM2STR402070N | AM2STR452070N | AM2STR502070N | AM2STR552070N | AM2STR602070N | AM2STR652070N |
| | 3 | | AM2STR403070N | AM2STR453070N | AM2STR503070N | AM2STR553070N | AM2STR603070N | AM2STR653070N |
| | 4 | | AM2STR403070N | AM2STR454070N | AM2STR504070N | AM2STR554070N | AM2STR604070N | AM2STR654070N |
| | 5 | | AM2STR405070N | AM2STR455070N | AM2STR505070N | AM2STR555070N | AM2STR605070N | AM2STR655070N |
| | 6 | | AM2STR406070N | AM2STR456070N | AM2STR506070N | AM2STR556070N | AM2STR606070N | AM2STR656070N |


STERIOSS IMPLANT Reference Code.

Angled Abutment (HEX)



- Used when a prosthetic's path adjustment is necessary at 15° or 25° axial angle
- Gold coloring for improved aesthetics
- 1.2mm Hex driver
- Recommended tightening torque : 30Ncm

Angled Abutment Guide




- 12 orientation choices with two types of hexes, H and V

| Cuff | | 4.0 | 4.5 | 5.0 | 5.5 | 6.0 | 6.5 |
|---------------|---|--------------|--------------|--------------|--------------|--------------|--------------|
| H TYPE | | | | | | | |
| 15° HEX | 1 | AM2AAR40115H | AM2AAR45115H | AM2AAR50115H | AM2AAR55115H | AM2AAR60115H | AM2AAR65115H |
| | 2 | AM2AAR40215H | AM2AAR45215H | AM2AAR50215H | AM2AAR55215H | AM2AAR60215H | AM2AAR65215H |
| | 3 | AM2AAR40315H | AM2AAR45315H | AM2AAR50315H | AM2AAR55315H | AM2AAR60315H | AM2AAR65315H |
| | 4 | AM2AAR40415H | AM2AAR45415H | AM2AAR50415H | AM2AAR55415H | AM2AAR60415H | AM2AAR65415H |
| | 5 | AM2AAR40515H | AM2AAR45515H | AM2AAR50515H | AM2AAR55515H | AM2AAR60515H | AM2AAR65515H |
| | 6 | AM2AAR40615H | AM2AAR45615H | AM2AAR50615H | AM2AAR55615H | AM2AAR60615H | AM2AAR65615H |
| 25° HEX | 1 | AM2AAR40125H | AM2AAR45125H | AM2AAR50125H | AM2AAR55125H | AM2AAR60125H | AM2AAR65125H |
| | 2 | AM2AAR40225H | AM2AAR45225H | AM2AAR50225H | AM2AAR55225H | AM2AAR60225H | AM2AAR65225H |
| | 3 | AM2AAR40325H | AM2AAR45325H | AM2AAR50325H | AM2AAR55325H | AM2AAR60325H | AM2AAR65325H |
| | 4 | AM2AAR40425H | AM2AAR45425H | AM2AAR50425H | AM2AAR55425H | AM2AAR60425H | AM2AAR65425H |
| | 5 | AM2AAR40525H | AM2AAR45525H | AM2AAR50525H | AM2AAR55525H | AM2AAR60525H | AM2AAR65525H |
| | 6 | AM2AAR40625H | AM2AAR45625H | AM2AAR50625H | AM2AAR55625H | AM2AAR60625H | AM2AAR65625H |
| V TYPE | | | | | | | |
| 15° HEX | 1 | AM2AAR40115V | AM2AAR45115V | AM2AAR45115V | AM2AAR45115V | AM2AAR60115V | AM2AAR65115V |
| | 2 | AM2AAR40215V | AM2AAR45215V | AM2AAR45215V | AM2AAR45215V | AM2AAR60215V | AM2AAR65215V |
| | 3 | AM2AAR40315V | AM2AAR45315V | AM2AAR45315V | AM2AAR45315V | AM2AAR60315V | AM2AAR65315V |
| | 4 | AM2AAR40415V | AM2AAR45415V | AM2AAR45415V | AM2AAR45415V | AM2AAR60415V | AM2AAR65415V |
| | 5 | AM2AAR40515V | AM2AAR45515V | AM2AAR45515V | AM2AAR45515V | AM2AAR60515V | AM2AAR65515V |
| | 6 | AM2AAR40615V | AM2AAR45615V | AM2AAR45615V | AM2AAR45615V | AM2AAR60615V | AM2AAR65615V |
| 25° HEX | 1 | AM2AAR40125V | AM2AAR45125V | AM2AAR45125V | AM2AAR45125V | AM2AAR45125V | AM2AAR45125V |
| | 2 | AM2AAR40225V | AM2AAR45225V | AM2AAR45225V | AM2AAR45225V | AM2AAR45225V | AM2AAR45225V |
| | 3 | AM2AAR40325V | AM2AAR45325V | AM2AAR45325V | AM2AAR45325V | AM2AAR45325V | AM2AAR45325V |
| | 4 | AM2AAR40425V | AM2AAR45425V | AM2AAR45425V | AM2AAR45425V | AM2AAR45425V | AM2AAR45425V |
| | 5 | AM2AAR40525V | AM2AAR45525V | AM2AAR45525V | AM2AAR45525V | AM2AAR45525V | AM2AAR45525V |
| | 6 | AM2AAR40625V | AM2AAR45625V | AM2AAR45625V | AM2AAR45625V | AM2AAR45625V | AM2AAR45625V |

STERIOSS IMPLANT Reference Code.


Angled Abutment (NON-HEX)



- Used when a prosthetic's path adjustment is necessary at 15° or 25° axial angle
- Gold coloring for improved aesthetics
- 1.2mm Hex driver
- Recommended tightening torque : 30Ncm

| | | Ø | 4.0 | 4.5 | 5.0 | 5.5 | 6.0 | 6.5 |
|-------------|------|---|--------------|--------------|--------------|--------------|--------------|--------------|
| 15° NON-HEX | Cuff | | | | | | | |
| | 1 | | AM2AAR40115N | AM2AAR45115N | AM2AAR50115N | AM2AAR55115N | AM2AAR60115N | AM2AAR65115N |
| | 2 | | AM2AAR40215N | AM2AAR45215N | AM2AAR50215N | AM2AAR55215N | AM2AAR60215N | AM2AAR65215N |
| | 3 | | AM2AAR40315N | AM2AAR45315N | AM2AAR50315N | AM2AAR55315N | AM2AAR60315N | AM2AAR65315N |
| | 4 | | AM2AAR40415N | AM2AAR45415N | AM2AAR50415N | AM2AAR55415N | AM2AAR60415N | AM2AAR65415N |
| | 5 | | AM2AAR40515N | AM2AAR45515N | AM2AAR50515N | AM2AAR55515N | AM2AAR60515N | AM2AAR65515N |
| 25° NON-HEX | 6 | | AM2AAR40615N | AM2AAR45615N | AM2AAR50615N | AM2AAR55615N | AM2AAR60615N | AM2AAR65615N |
| | 1 | | AM2AAR40125N | AM2AAR45125N | AM2AAR50125N | AM2AAR55125N | AM2AAR60125N | AM2AAR65125N |
| | 2 | | AM2AAR40225N | AM2AAR45225N | AM2AAR50225N | AM2AAR55225N | AM2AAR60225N | AM2AAR65225N |
| | 3 | | AM2AAR40325N | AM2AAR45325N | AM2AAR50325N | AM2AAR55325N | AM2AAR60325N | AM2AAR65325N |
| | 4 | | AM2AAR40425N | AM2AAR45425N | AM2AAR50425N | AM2AAR55425N | AM2AAR60425N | AM2AAR65425N |
| | 5 | | AM2AAR40525N | AM2AAR45525N | AM2AAR50525N | AM2AAR55525N | AM2AAR60525N | AM2AAR65525N |
| | 6 | | AM2AAR40625N | AM2AAR45625N | AM2AAR50625N | AM2AAR55625N | AM2AAR60625N | AM2AAR65625N |

Milling Abutment



| | | Ø | 4.0 | 4.5 | 5.0 | 6.0 |
|---------|------|---|-------------|-------------|-------------|-------------|
| HEX | Cuff | | | | | |
| | 1.5 | | AM2MAR4015H | AM2MAR4515H | AM2MAR5015H | AM2MAR6015H |
| NON-HEX | 3 | | AM2MAR4030H | AM2MAR4530H | AM2MAR5030H | AM2MAR6030H |
| | 1.5 | | AM2MAR4015N | AM2MAR4515N | AM2MAR5015N | AM2MAR6015N |
| | 3 | | AM2MAR4030N | AM2MAR4530N | AM2MAR5030N | AM2MAR6030N |

- Used when an abutment's path must be altered or a prosthetic's margin area must be customized
- Gold coloring for improved aesthetics
- 1.2mm Hex driver
- Recommended tightening torque : 30Ncm

STERIOSS IMPLANT Reference Code.

CCM Abutment



HEX

ECCRH401C

NON-HEX


ECCRN401C

Abutment Screw



AM2ASR

Temporary Abutment



- Used in producing temporary prosthetics
- Structure enabling easy customization
- 1.2mm Hex driver
- Recommended tightening torque : 20Ncm

Type Ø

4.5

HEX


AM2TAR 4510H

NON-HEX

AM2TAR4510N

STERIOSS IMPLANT Reference Code.


Pick-up Impression Coping



| | | 4.0 | 4.5 | 5.0 | 5.5 | 6.0 | 6.5 |
|---------|-------|------------|------------|------------|------------|------------|------------|
| HEX | SHORT | AM2PIR40SH | AM2PIR45SH | AM2PIR50SH | AM2PIR55SH | AM2PIR60SH | AM2PIR65SH |
| | LONG | AM2PIR40LH | AM2PIR45LH | AM2PIR50LH | AM2PIR55LH | AM2PIR60LH | AM2PIR65LH |
| NON-HEX | SHORT | AM2PIR40SN | AM2PIR45SN | AM2PIR50SN | AM2PIR55SN | AM2PIR60SN | AM2PIR65SN |
| | LONG | AM2PIR40LN | AM2PIR45LN | AM2PIR50LN | AM2PIR55LN | AM2PIR60LN | AM2PIR65LN |


- Takes impression using open tray technique
- Superior impression stability with an extended flat & groove
- 1.2mm Hex driver

Impression Coping Guide



STERIOSS IMPLANT Reference Code.

Transfer Impression Coping



- Takes impression using closed tray technique
- 1.2mm Hex driver

| | | Ø | 4.0 | 4.5 | 5.0 | 5.5 | 6.0 | 6.5 |
|---------|-------|---|------------|------------|------------|------------|------------|------------|
| HEX | SHORT | | AM2TIR40SH | AM2TIR45SH | AM2TIR50SH | AM2TIR55SH | AM2TIR60SH | AM2TIR65SH |
| | LONG | | AM2TIR40LH | AM2TIR45LH | AM2TIR50LH | AM2TIR55LH | AM2TIR60LH | AM2TIR65LH |
| NON-HEX | SHORT | | AM2TIR40SN | AM2TIR45SN | AM2TIR50SN | AM2TIR55SN | AM2TIR60SN | AM2TIR65SN |
| | LONG | | AM2TIR40LN | AM2TIR45LN | AM2TIR50LN | AM2TIR55LN | AM2TIR60LN | AM2TIR65LN |

Fixture Lab Analog



- Achieves fixture in oral cavity on a working model

| Ø | 3.5 | 4.0 |
|---|--------|--------|
| | FLAR35 | FLAR40 |


Abutment Lab Analog



| Height | Ø | 4.0 | 4.5 | 5.0 | 6.0 |
|--------|---|---------|---------|---------|---------|
| 4 | | ALA4040 | ALA4440 | ALA5040 | ALA6040 |
| 5.5 | | ALA4055 | ALA4455 | ALA5055 | ALA6055 |
| 7 | | ALA4070 | ALA4470 | ALA5070 | ALA6070 |

PROSTHETICS FLOW DIAGRAMS.

Screw & Cement Retained Restoration



STERIOSS IMPLANT Reference Code.


Esthetic Solid Abutment



- Used in screw type prosthetics in multiple cases
- Gold coloring for improved aesthetics
- Recommended tightening torque : 30Ncm

| Cuff | Ø | 4.8 |
|------|---|----------|
| 1 | | AM2ESR10 |
| 2 | | AM2ESR20 |
| 3 | | AM2ESR30 |
| 4 | | AM2ESR40 |
| 5 | | AM2ESR50 |
| 6 | | AM2ESR60 |

Esthetic Angled Abutment



- Used in implant path compensation in screw retained multiple case
- Up to 60° path compensation (two implant standard)
- Gold coloring for improved aesthetics
- Recommended tightening torque : 30Ncm
- 1.2mm Hex driver

| Cuff | Ø | 4.8 |
|------|---|------------|
| 2 | | AM2EAR2017 |
| 3 | | AM2EAR3017 |
| 4 | | AM2EAR4017 |
| 17° | | |
| 3 | | AM2EAR3030 |
| 4 | | AM2EAR4030 |
| 5 | | AM2EAR5030 |
| 30° | | |

STERIOSS IMPLANT Reference Code.

Esthetic Healing Cap




- Used when protecting a convertible abutment in the oral cavity and minimizing foreign body sensation for the patient
- Recommended tightening torque : 20Ncm
- 1.2mm Hex driver

Esthetic Pick-up Impression Coping



- Takes impression using open tray
- Superior impression stability with an extended flat & groove
- 1.2mm Hex driver

Esthetic Temporary Abutment



- Used in producing temporary prosthetics
- Structure enabling easy customization and minimizing indication restrictions
- 1.2mm Hex driver
- Recommended tightening torque : 20Ncm


STERIOSS IMPLANT Reference Code.

Esthetic Plastic Abutment



- Used in producing screw maintenance prosthetics
- Produces prosthetics after casting with dentalgrade alloy (gold, non-precious metals) after customization
- 1.2mm Hex driver
- Recommended tightening torque : 20Ncm


Esthetic Lab Analog



- Achieves Esthetic Solid Abutment of the oral cavity on a working model

PROSTHETICS FLOW DIAGRAMS.

Overdenture Restoration



STERIOSS IMPLANT Reference Code.

Ball Abutment



- Used in creating stud type overdenture prosthetics
- Gold coloring for improved aesthetics
- Rotator Abutment driver
- Recommended tightening torque : 30Ncm

| Cuff | Ø | 4.0 |
|------|---|----------|
| 0 | | AM2BAR00 |
| 1 | | AM2BAR01 |
| 2 | | AM2BAR02 |
| 3 | | AM2BAR03 |
| 4 | | AM2BAR04 |
| 5 | | AM2BAR05 |
| 6 | | AM2BAR06 |

Ball Cap



- Used in creating stud type overdenture prosthetics

| Ø | 4.0 |
|---|-------|
| | AM2BC |


Ball Lab Analog



| Ø | 3.5 |
|---|-------|
| | BLA35 |

PROSTHETICS FLOW DIAGRAMS.

Overdenture Restoration



STERIOSS IMPLANT Reference Code.

Rotator Abutment



- Used in creating Rotator type overdenture prosthetics
- Gold coloring for improved aesthetics
- Rotator Abutment driver
- Recommended tightening torque : 30Ncm

| Ø | |
|------|-----------|
| Cuff | |
| 1 | LOGSR3710 |
| 2 | LOGSR3720 |
| 3 | LOGSR3730 |
| 4 | LOGSR3740 |
| 5 | LOGSR3750 |
| 6 | LOGSR3760 |
| 7 | LOGSR3770 |
| 8 | LOGSR3780 |

Impression Cap



- Used in creating Rotator type overdenture prosthetics

| Ø | |
|-----|-------|
| 4.0 | ROTIP |

Rotator Analog



| Ø | |
|-----|----------|
| 4.0 | ROTAL 40 |

- Achieves Rotator Abutment of the oral cavity on a working model


STERIOSS IMPLANT TL(Tissue Level) SYSTEM Reference Code.

FIXTURE starts from 6mm and has a maximum of 13mm every 1mm.

3.5 / 4.0 / 4.5 / 5.0




5.0 / 6.0



FISR


Fixture System Regular(FISR)
5.0 Prosthetic Connection



Fixture System Regular(FISR)
5.0 Prosthetic Connection is
compatible with all Straumann
ITI Regular Neck 4.8mm Implants.

FISW

Fixture System Wide(FISW)
6.0 Prosthetic Connection



Fixture System Wide(FISW)
6.5 Prosthetic Connection is
compatible with all Straumann
ITI Wide Neck 6.5mm Implants.


TL(Tissue Level) SYSTEM Reference Code.

FIXTURE starts from 6mm and has a maximum of 13mm every 1mm.



FISR

| Cuff Ø | 3.5 | 4.0 | 4.5 | 5.0 |
|--------|--------------|--------------|--------------|--------------|
| 3 | | | | FISR 513030 |
| 4 | | FISR 423040 | FISR 463040 | FISR 513040 |
| 5 | | FISR 423050 | FISR 463050 | FISR 513050 |
| 6 | | FISR 423060 | FISR 463060 | FISR 513060 |
| 7 | FISR 383070 | FISR 423070 | FISR 463070 | FISR 513070 |
| 8 | FISR 383080 | FISR 423080 | FISR 463080 | FISR 513080 |
| 9 | FISR 383090 | FISR 423090 | FISR 463090 | FISR 513090 |
| 10 | FISR 3830100 | FISR 4230100 | FISR 4630100 | FISR 5130100 |
| 11 | FISR 3830110 | FISR 4230110 | FISR 4630110 | FISR 5130110 |
| 12 | FISR 3830120 | FISR 4230120 | FISR 4630120 | FISR 5130120 |
| 13 | FISR 3830130 | FISR 4230130 | FISR 4630130 | FISR 5130130 |




FISW

| Cuff Ø | 5.0 | 6.0 |
|--------|--------------|--------------|
| 3 | FISW 513030 | FISW 603030 |
| 4 | FISW 513040 | FISW 603040 |
| 5 | FISW 513050 | FISW 603050 |
| 6 | FISW 513060 | FISW 603060 |
| 7 | FISW 513070 | FISW 603070 |
| 8 | FISW 513080 | FISW 603080 |
| 9 | FISW 513090 | FISW 603090 |
| 10 | FISW 5130100 | FISW 6030100 |
| 11 | FISW 5130110 | FISW 6030110 |
| 12 | FISW 5130120 | FISW 6030120 |
| 13 | FISW 5130130 | FISW 6030130 |

Dia 3.5 / 4.0 / 4.5 / 5.0



Dia 5.0 / 6.0



STERIOSS IMPLANT Reference Code.

Connection Guide

Cover Screw

Headless
Cover Screw



SIHLS

Cover Screw




SICSR 48

WD Cover Screw



SICSW 65

Healing Abutment




- 1.2mm Hex driver
- Recommended tightening torque : 5~8Ncm

| Height | Ø | SIHARDP |
|--------|---|----------|
| 4.8 | | SIHARDP1 |
| | | SIHARDP2 |
| | | SIHARDP3 |
| | | SIHARDP4 |

| Height | Ø | SIHAWDP |
|--------|---|----------|
| 6.5 | | SIHAWDP1 |
| | | SIHAWDP2 |
| | | SIHAWDP3 |
| | | SIHAWDP4 |


PROSTHETICS FLOW DIAGRAMS.

Screw & Cement Retained Restoration



STERIOSS IMPLANT Reference Code.

Straight Abutment



Cuff / Ø 4.8 RP


| | |
|---|-------------|
| 0 | SISTARDP060 |
| 1 | SISTARDP160 |
| 2 | SISTARDP260 |
| 3 | SISTARDP360 |
| 4 | SISTARDP460 |

Cuff / Ø 6.5 WP

| | |
|---|-------------|
| 0 | SISTAWDP060 |
| 1 | SISTAWDP160 |
| 2 | SISTAWDP260 |
| 3 | SISTAWDP360 |
| 4 | SISTAWDP460 |

STERIOSS IMPLANT Reference Code.

Angled Abutment



| | 15° | 25° |
|------|----------------|----------------|
| (RP) | SIAAR DP 06150 | SIAAR DP 06250 |
| | 16150 | 16250 |
| | 26150 | 26250 |
| | 36150 | 36250 |
| | 46150 | 46250 |

| | 15° | 25° |
|------|----------------|----------------|
| (WP) | SIAAW DP 06150 | SIAAW DP 06250 |
| | 16150 | 16250 |
| | 26150 | 26250 |
| | 36150 | 36250 |
| | 46150 | 46250 |


CCM Abutment



CCM Abutment

| | OCTA | NON OCTA |
|------|------------|------------|
| (RP) | ICCORH481D | ICCORN481D |
| (WP) | ITW-O | ITW-N |

Temporary Abutment



RP

SITARDPO

WP

SITAWDPO

- Used in producing temporary prosthetics
- Structure enabling easy customization
- 1.2mm Hex driver
- Recommended tightening torque : 20Ncm

Pick-up Impression Coping

RP


SIPIRDPSO

SIPIRDPLO


WP

SIPIWDPSO

SIPIWDPLO



Transfer Impression Coping



- Takes impression using closed tray technique
- 1.2mm Hex driver

RP

SITIRDPSO

SITIRDPLO

WP

SITIWDPSO

SITIWDPLO

STERIOSS IMPLANT Reference Code.

Solid Abutment



Hight

RP

WP

3mm

SISARDP3

SISAWDP3

4mm

SISARDP4

SISAWDP4

5mm

SISARDP5

SISAWDP5

6mm

SISARDP6

SISAWDP6

7mm

SISARDP7

SISAWDP7

Plastic Imperssion Coping



RP

WP

SIPCRDP7

SIPCWDP7

Fixture Lab Analog



RP

WP

SIFLARDP

SIFLAWDP

- Achieves fixture in oral cavity on a working model

Solid Lab Analog



RP

WP

SIFALRDP3

SIFALWDP3

SIFALRDP4

SIFALWDP4

SIFALRDP5

SIFALWDP5

SIFALRDP6


SIFALWDP6

SIFALRDP7

SIFALWDP7


PROSTHETICS FLOW DIAGRAMS.

Screw & Cement Retained Restoration




PROSTHETICS FLOW DIAGRAMS.

Overdenture & Cement Retained Restoration



STERIOSS IMPLANT Reference Code.

Rotator Abutment



| Ø | |
|------|------------|
| Cuff | |
| 0 | LGSSSR4800 |
| 1 | LGSSSR4810 |
| 2 | LGSSSR4820 |
| 3 | LGSSSR4830 |
| 4 | LGSSSR4840 |

- Used in creating stud type overdenture prosthetics
- Gold coloring for improved aesthetics
- Ball Abutment driver
- Recommended tightening torque : 30Ncm

Impression Cap



| | |
|---|-------|
| Ø | 4.0 |
| | ROTIP |

- Used in creating Rotator type overdenture prosthetics

Rotator Lab Analog



| | |
|---|----------|
| Ø | 4.0 |
| | ROTAL 40 |

- Achieves Rotator Abutment of the oral cavity on a working model

Kit Rendering




Surgical Kit Lay-out

Point Drill



TPD16

Parallel Pin




TTP

Drill Extension



TDE



Lindermann Drill



TLD

Crest Burr



TCB

BL H.P Connector



THC 25L



THC 25S

BL Torque Connector



TWTC 25L



TWTC 25M

Hand Driver



THD12L



THD12S

Spread Drill

TSPD35



TSPD40



TSPD45



TSPD50



TSPD60



Torque Wrench




STW01


Surgical Kit

Step Drill


Ø2.2




F3.5 (Ø3.5)




F4.0 (Ø4.0)




F4.5 (Ø4.4)



F5.0 (Ø4.9)



F6.0 (Ø5.8)




Full Surgical Kit STCK1

| Point Drill [TPD16] | Lindermann Drill [TLD] | Crest Bur [TCB] | Drill Extension [TDE] |
|-----------------------|--|--|--|
| Step Drill | [TSD2206] [TSD2207] [TSD2208] [TSD2209] [TSD2210] [TSD2211] [TSD2212] [TSD2213] | [TSD3506] [TSD3507] [TSD3508] [TSD3509] [TSD3510] [TSD3511] [TSD3512] [TSD3513] | [TSD4006] [TSD4007] [TSD4008] [TSD4009] [TSD4010] [TSD4011] [TSD4012] [TSD4013] |
| Spreader Drill | [TSPD35] [TSPD40] [TSPD45] [TSPD50] [TSPD60] | Parallel Pin H.P Connector Ratchet Connector | 3EA Long THC25L Short THC25S Long TWTC25L Short TWTC25M |
| Torque Wrench [STW01] | | Hand Driver | Long THD12L Short THD12S |

Standard Surgical Kit STCK2

| Point Drill [TPD16] | Lindermann Drill [TLD] | Crest Bur [TCB] | Drill Extension [TDE] |
|---------------------|---|---|---|
| Step Drill | [TSD2206] [TSD2207] [TSD2208] [TSD2209] [TSD2210] | [TSD3506] [TSD3507] [TSD3508] [TSD3509] [TSD3510] | [TSD4006] [TSD4007] [TSD4008] [TSD4009] [TSD4010] |
| Parallel Pin | 3EA | Torque Wrench [STW01] | |
| H.P Connector | Long THC25L Short THC25S | | |
| Ratchet Connector | Long TWTC25L Short TWTC25M | | |
| Hand Driver | Long THD12L Short THD12S | | |

Full Surgical Kit STCK3



Point Drill [TPD16] **Lindermann Drill** [TLD] **Crest Bur** [TCB] **Drill Extension** [TDE]

| | | | | | | |
|-------------------|-----------|-----------|-----------|-----------|-----------|-----------|
| Step Drill | [TSD2203] | [TSD3503] | [TSD4003] | [TSD4503] | [TSD5003] | [TSD6003] |
| | [TSD2204] | [TSD3504] | [TSD4004] | [TSD4504] | [TSD5004] | [TSD6004] |
| | [TSD2205] | [TSD3505] | [TSD4005] | [TSD4505] | [TSD5005] | [TSD6005] |
| | [TSD2206] | [TSD3506] | [TSD4006] | [TSD4506] | [TSD5006] | [TSD6006] |
| | [TSD2207] | [TSD3507] | [TSD4007] | [TSD4507] | [TSD5007] | [TSD6007] |
| | [TSD2208] | [TSD3508] | [TSD4008] | [TSD4508] | [TSD5008] | [TSD6008] |
| | [TSD2209] | [TSD3509] | [TSD4009] | [TSD4509] | [TSD5009] | [TSD6009] |
| | [TSD2210] | [TSD3510] | [TSD4010] | [TSD4510] | [TSD5010] | [TSD6010] |

Spreader Drill [TSPD35]
[TSPD40]
[TSPD45]
[TSPD50]
[TSPD60]

Parallel Pin 3EA


H.P. Connector Long THC25L
Short THC25S

Ratchet Connector Long TWTC25L
Short TWTC25M

Torque Wrench [STW01]

Hand Driver Long THD12L
Short THD12S

Short Drill Kit SDK01



Implant Torque Wrench ITW01



Sterioss Multi Prosthetic Kit **MPK01**



1.2Hex Machine driver

| | |
|-------|-----|
| Short | 1ea |
| Long | 1ea |

1.2Hex Hand driver

| | |
|-------|-----|
| Short | 1ea |
| Long | 1ea |

1.27Hex Hand driver

| | |
|-------|-----|
| Short | 1ea |
| Long | 1ea |

Strauman Hand driver

| | |
|-------|-----|
| Short | 1ea |
| Long | 1ea |

Nobel Hand driver

| | |
|-------|-----|
| Short | 1ea |
| Long | 1ea |

Slot driver

| | |
|------|-----|
| Long | 1ea |
|------|-----|

Connector


1ea

Torque Wrench


1ea

Drilling Sequence


F3.5 (Ø3.8)



F3.5 Spread Drill




F4.0 (Ø4.2)




F4.0 Spread Drill



F4.5 (Ø4.6)





F4.5 Spread Drill




Drilling Sequence

F5.0 (Ø5.1)



F6.0 (Ø6.0)



High Power System Fourth Generation New Implant Engine

LED Light

Powerful optic motor
torque function
Touch control LCD screen
User friendly foot control

Model : STERIOSS PRO
Intuitive powerful safe and reliable

Maxim Power 80Ncm

Technical:

- Maximun power : 150W
- Power supply : AC230V 50/60Hz
- Motor speed : 300--40000rpm
- Adustable torque : 5Ncm--80Ncm
- Water supply volume :
4 blocks adjustable,
150ml/min
- Light source illumination : 30000lux
- Motor : surgical brushless motor
- Volume(cm) : 42x38x26cm
- Weight(kg)/PCS : 5.5kg
- Packing material : Aluminum box
- QTY per master carton : 1 pc/ctn

- Full touch big screen

The large size screen ensures that the treatment parameters are visible at any time

- Graphical planting process selection, a variety of parameters at a glance

Every planting process is set up to solve the thdious setting

- ad hoc separate irrigation procedures for ease of use

Cooling and flushing of inner and outer channel

- Built in water supply system with small size, good performance and convenient convection

4 shifts water supply up to 150ml/min

- a new generation of motors with torque of up to 80N.cm

High performance brushless motor, powerful touque output, highly efficient chip structure

- motors and cables can be sterilized at high temperature and pressure

Lighter and shorter motor, long work no fatigue

- LED lighting helps clinicians focus More on the treatment area and make surgery more accurate

High brightness LED light source with brightness up to 30000lux


- illumination and sound can be adjusted

Give doctors more autonomy in setting up options

- Multifunctional pedal controller

Mseptic operation can be achieved, giving the maximum freedom of the handle

The handle can be easily connected to the pedal controller



User Manual(Fixture)

1) Product description

Sterioss Implant is made of titanium and grafted in the maxillary or mandibular alveolar bone for surgical treatment. It plays the role of natural dental root. Refer to the manual or the catalogue or our website (www.steriossimplant.com) for detail. See the product label for the product code, specifications, manufacturing date, and expiration date.

2) Use

Should be used in accordance with the general implant guide.

3) Indications

- Generally good health, healthy oral cavity
- Predicted normal healing of wounds
- Jaw bone that has stopped growing
- Good oral cavity hygiene
- Available bone sufficient for vertical, mesio/distal, and buccal/lingual space
- Patient's full recognition of the risk in the course of performing each implant treatment

4) Contraindications

Following medical conditions may occur failure in cases when:

- Insufficient bones or poor bone quality to achieve implant stability
- Poor oral hygiene or pathological infection to alveolar bones
- Patient with uncontrollable diabetes, tissue disease influencing bone or wound treatment; heavy smoker or alcoholic
- Heavy smoking, tobacco abuse
- Any patient who is not suitable for operation

5) Warnings

The selection of inappropriate patients and operation methods can cause implant failure or loss of bone supporting the implant. Sterioss implants must not be used for purposes other than the recommended use and must not be remodeled. Implant mobility, bone loss, and chronic infection can result in failure of the implant surgery.

6) Precautions

The surgical technology of dental implant involves an expert, complex procedure. Formal training is required to perform implant surgery. Determine the local anatomy and suitability of the available bone for implant placement. Prepare the implant considering the expected situations and cautions. Visual inspections as well as panoramic and periapical radiographs are essential to determine anatomical landmarks, occlusal conditions, periodontal status, and the adequacy of the bone. Adequate radiographs, direct palpation, and visual inspection of the implant site are necessary prior to treatment, planning and use of Sterioss implant.

7) Adverse reactions

A few problems may occur after the operation (loss of implant stability, damage of prosthesis, etc.) deficient quality and quantity of the remaining bone, infection, inferior oral hygiene or uncooperativeness of patient, implant mobility, partial deterioration of tissue, and improper position or arrangement of implants may cause the above mentioned problems.

8) Surgical complications

The implant procedure has risk, including localized swelling, dehiscence, tenderness of short duration, edema, hematoma, or bleeding. Numbness of the lower lip and chin region following lower jaw surgery, and of the tissue beside the nose following upper jaw surgery, is a possible side-effect of the surgery. Though it would most probably be of a temporary nature, in very rare cases, the numbness has been permanent. Gingival-mucosal (gum tissue) ulceration, tissue reaction, or infection may occur, but generally responds to local care.

9) Sterilization and handling

The implant has been cleaned and sterilized by gamma irradiation and is ready to use. The sterilized product must be used in a sterilized environment with sterilized tools. If the packaging is damaged, or if the expiration date has passed, do not use the product. Expired or contaminated product must not be re-sterilized; they must be disposed of.

10) Storage Condition

Keep the product in a dry place at room temperature (1°C~20°C). Keep away from direct sunlight. This product is valid for 3 years from the date of manufacture, please observe the validity date.

11) Cautions for use

Determine the number, angles, and positions of the implants to be grafted considering a minimum gap of 3mm between implant. Local anesthesia is used; an anesthetic solution should be sufficient to anesthetize the periosteum surface. Conventional implant treatment should be applied. When using engine driver, the recommended speed is 15rpm, and maximum torque, 35Ncm. Once the maximum torque is reached, detach the driver mounted in the hand piece. Afterward, connect the torque wrench and perform grafting while rotating the hand-driver clockwise. It must be checked the information on length, diameter and depth for placement before use.

User Manual(Abutment)

1) Product description

Sterioss Implant is a dental implant materials and consist of the Titanium. Refer to manual, catalog or our website (www.zerosimplant.com) for details. For the product code, specification, manufacturing date, and expiration date see the product label.

2) Use

Should be used in accordance with the general implant guide.

3) General Precaution

The surgical technology of dental implantation requires an expert, it is a complex procedure which requires formal training to perform implantation.

4) Important

It is important to look at the anatomy and suitability of the available bone for implant placement. Prepare the implant considering the expected situations and cautions. Visual inspection as well as panoramic and periapical radiographs are essential to determine anatomical landmarks, occlusal conditions, periodontal status, and adequacy of bone. Lateral cephalometric radiographs, CT scans, and tomograms may also be beneficial. In particular, the exact implant fixture can be assembled to prepare the abutment and the prosthetic components.

5) Indication

- Partial or full edentulous patients
- Periodontitis
- Dental caries
- Accidental tooth damage
- Congenital tooth defects
- Tooth loss due to other maxillofacial diseases

6) Contraindications

Contraindications include following, but are not limited to:

- Hemophilia Patient
- Patient experiencing difficulties related to bone and wound treatment
- Patient with uncontrollable diabetes, tissue disease influencing bone or wound treatment; heavy smoker or alcoholic
- Patient whose immunity system is inactivated due to chemical therapy and radiation therapy
- Patient with oral infection or inflammation (Improper oral hygiene, bruxism)
- Patient with untreatable occlusion/joint disorder, insufficient dental arch space
- Any patient who is not suitable for operation

7) Procedural Precautions (Surgery)

The Implant operation requires high accuracy and careful attention, we must try to minimize damage to the cell tissue and pay special attention to the temperature, surgical trauma, and/or removal of the source of contamination and infection.

8) Procedural Precautions (Prosthetics)

The prosthetic structure is small; make sure it is neither swallowed nor inhaled by the patient. Angled abutments are not recommended for placement in the posterior region of the mouth due to limitations of implant strength. Stress distribution is especially important in implant operation as well as the fit of prosthesis and abutment on bridge, also the occlusal stability. Avoid using excessive force horizontally especially during immediate implantation. For the prosthesis whose substructure is made of gold alloy, gold should be used appropriately. Operate prosthesis after enough healing period.

9) Cautions for Patients

Do not apply excessive stress on the replacement until the last prosthesis is placed.

10) Side Effect

Side Effect

These Problems may occur after implantation (loss of implant stability, loss of prosthesis, etc). Deficient quality and quantity of remaining bone, infection, inferior oral hygiene or uncooperativeness of patient, implant mobility, partial deterioration of tissue, and improper position and arrangement of implants can cause instability.

11) Warning

Read carefully instruction for use before use.

Using Sterioss Implant safely and effectively requires special training, since the surgical techniques involved in the dental implant operation are highly specialized and very complex. The selection of inappropriate patients and operation methods can cause implant failures or loss of bone supporting the implant. Sterioss Implant must not be used for purposes other than the recommended use. Dental implants must never be remodeled. If the implant is contaminated by the patient's bodily fluids, it cannot be used in other patients. The ceramic abutment needs a special manufacturing process; the technician should be specifically trained in this process.

12) Sterility

The unsterilized prosthetic components must be sterilized in an autoclave at 132°C for 15 minutes before use. After the steam sterilization, the abutments should be dried for 15 minutes before use. This product is disposable medical device product in any case should not reuse.

13) Storage Conditions

Keep in cool(1°C~30°C) and dry place. Keep Out of direct sunlight This product is valid for semi-permanent before open.

Product compatibility list

GBR Product

Allograft

PREMIUM



● Renew Oss The Black (syr.)

Power: 0.2~1mm
Volume : 0.3cc / 1cc



● Do Bone (syr.)

Power: S:0.25~0.71mm /
L:0.71~1.6mm
Volume : 0.3cc / 1cc

HIT



● Renew Oss (Syr.)

Power: 0.2~1mm
Volume : 0.3cc / 0.5cc / 1cc

NEW



● Titan Bone (syr.)

Volume : 0.3cc / 0.6cc / 1.1cc



● Allobone Plus

Volume : 0.25cc / 0.5cc / 1cc



● Oragraft

Volume : 0.25cc / 0.5cc / 1.2cc

DBM/BMP

<DBM>



● Orthoblast II (Gel)

Volume : 0.5cc / 5cc / 3cc

<DBM>



● DynaGraft II

Volume : 0.5cc / 1cc / 2.5cc / 5cc

HIT



● Demios

Volume : 0.25cc / 0.5cc / 1cc

<DBM>

PREMIUM

<BMP>



● Osteosparx-C

Volume : 1cc / 3cc / 5cc / 10cc

PREMIUM

<BMP>



● Allomatrix

Volume : 1cc

<BMP>



● Rafugen BMP2

Volume : 0.1cc / 0.3c / 0.5cc

Synthetic

HIT



● Excelos Inject

Volume : 0.45g / 0.75g / 1.5g



● Sorbone Syringe

Power: M:0.5~1mm / L:1~2mm
Volume : 0.25cc / 0.5cc / 1cc

HIT



● Sorbone

Power: S:0.3~0.5mm / M:0.5~1mm / L:1~2mm
Volume : 0.25g / 0.5g / 2g



● SynCera II

Power: S:0.4~1mm / L:1~2mm
Volume : 0.25g / 0.5g / 1g

NEW



● Bongros

Power: M:0.6~1mm / L:1~3mm
Volume : 0.25g(0.4cc) / 0.5g(0.8cc)



● Straumann Synthetic

Power: M:0.6~1mm / L:1~3mm
Volume : 0.25g / 0.5g / 1g

Xenograft



● NeoOss

Volume : 0.25g (≒0.5cc)



● NeoOss Plug (Porcine+Collagen)

Size: Ø6~10mm / Ø10x10mm

NEW



● Titan-X

Power: S:0.4~1.0mm
Volume : 0.1g / 0.25g / 1g / 2g



● Bio-Oss

Volume : S:0.25~1mm
(0.25g / 0.5g / 1g / 2g)
L:1~2mm(0.5g / 1g / 2g)



● Bio-Oss Pen

Volume : S:0.25~1mm(0.5cc / 1cc)
L:1~2mm(15cc)



● Bio-Oss Collagen

Volume : 100mg / 250mg



● Straumann Xeno

Power : S(0.2~1mm) L:1.0~2.0mm
Volume : 0.15g / 0.25g / 0.5g / 1g / 2g



● InduCera

Power : S(0.25~1mm) / L:1~2mm
Volume : 0.15g / 0.25g / 0.5g / 1g



● Cerabone

Power : S:0.5~1mm / L:1~2mm
Volume : 0.5cc / 1cc / 2cc / 5cc

Hemostasis / Absorbable



● Atelo Plug

Size : S:8x25mm(5ea) /
M:15x25mm(5ea)



● AteloCare (Sheet)

Size : 25x75x0.8mm(5ea)

HIT



● CollaDerm (Plug)

Size : S:8x25mm(5ea)
M:15x25mm(5ea)



● Colla Tape

Size : 25x75mm(10ea)



● Colla Plug

Size : 10x20mm(10ea)



● Colla Cote

Size : 10x20mm(10ea)

Membrane

NEW



■ Diaderm M

Size : 15x30mm



● Redura

Size : 15x20mm / 20x30mm



● Lyopalant

Size : 15x30mm



● CollaGuide

Size : 15x30mm / 15x30mm / 30x40mm



● LysoGide

Size : 10x20mm / 20x30mm



● Membra-Gen

Size : 15x18mm / 15x25mm

GBR Product



● Bio-Gide

Size : 13x25mm / 25x25mm
30x40mm



● CGDerm

Size : 10x20mm / 10x40mm
20x20mm / 20x40mm
Thickness : 10/20/30



● Remaix

Size : 15x20mm / 25x30mm
30x40mm



● Bioflexgide

Volume/Size : 10x20mm



● Jason Membrane

Volume/Size : 15x20mm / 20x30mm
30x40mm

Membrane / Non Absorbable



● Cytoflex TEF-Guard

Size : 12x24mm



● Cytoflex Mesh

Size : 25x30mm



● Renew Mesh-Ti

Size : 15x22mm / 22x36mm / 44x44mm

TITANIUM-REINFORCED MICRO PORE PTFE MEMBRANE



Instruments System

001

Total Haenaem Bur Set



002

Haenaem Bur Kit for Sinus



003

Haenaem Bur Kit for Expander



004

Bone Expander Kit



005

Double A Guide Kit



006

Total Sinus Kit



007

Crestal Approach Sinus Kit



008

One Drilling System Kit



009

Total Remove Kit



010

Solid Screw Kit for GBR



011

Bone Collector Kit



012

V-Bone Collector Kit



013

Initial Kit



014

Bone Mill Kit



015

Tissue Punch Kit



016

Round Bur Kit



017

Trephine Kit



#001

Total HaeNaem Bur Set





World Patent & Design by HaeNaem “ZERO” Bone Loss Drill

Sinus Auto Grafting/Ridge Expansion/ D4->D2 Bone densification at once with simple drilling



1. Maxillary Sinus Autografting is secure only by Drilling
2. Improves bone density from D4 to D2 which is weak due to drilling
3. Excellent Ridge expansion effect only by drilling
4. Improves bone density from D4 to D2 which is weak bone density by drilling.
5. Safe surgical operation is possible due to the clockwise drilling like the normal drill.
6. Easier and safer operation using the stoppers.



Comparative experiment

1) Drilling test in D3 Bone block



Normal Drill


The bone particles are discharged to the opposite side way (Back side).



Haenaem Bur Drill

The bone particles gather in the front direction and densification occurs.

2) Liquid experiment





Normal Drill







Haenaem Bur Drill

3) Normal Drill / Haenaem Bur Drill



Normal: Side perforation
HN Bur: Expansion without side perforation.

Easy & Simple Operation




Sinus

• PILOT DRILL



• SINUS DRILL (600~1200 RPM, 600 RPM when proceed perforation of catilage bone)



- A structure in which 70% of the force is transmitted in the traveling direction and 30% of the force is transmitted in the lateral direction.
- Sinus Autografting / Bone quality improvement possible only by drilling.
- During drilling, the remaining bones and cell lines rise to the maxilla at the same time as the maxillary sinus and lower cartilage are perforated.

<Clinical Data>

1) No.16 Sinus Lift



2) No.15, 16 Sinus Lift immediately after extraction



3) No.14, 15 Sinus Lift immediately after extraction



4) No.6 Sinus Lift immediately after extraction



5) No.4, 5, 6 Sinus Lift



Expander

• PILOT DRILL



• EXPANDER DRILL (600~1200 RPM, Depending on bone density)




- A structure in which 20% of the force is transmitted in the traveling direction and 80% of the force is transmitted in the lateral direct
- Excellent Ridge Expansion effect in narrow alveolar bone, extraction and septum only by drilling.
→Improved bone density from weak D4 to D2 bone.
- The drill does not slip during drilling, so the foresight of surgery is high.

<Clinical Data>


1) Expander




2) Septum Expander Case 1




3) Septum Expander Case 2



4) Septum Expander Case 3



5) Septum Expander Case 4







Residual bone 5mm+4.0 Drilling sequence for implant placement




※ All drilling should be based on pumping movements that repeat Up and Down and sense of pushing bones.

Drilling sequence by fixture size


| | |
|------------------------------|---|
| Drilling RPM | 600~1200 RPM (600 RPM when proceed perforation of cartilage bone) |
| Using Artificial Bone | Water OFF with Final Drill, RPM 50~100 |
| Normal Bone |  Drilling sequence |
| Soft Bone |   Drilling sequence |
| Hard Bone |  After drilling to half of the next optional drill, place the fixture |

[1] Sinus sequence


1) 4.0 fixture Placement




2) 4.5 fixture Placement



3) 5.0 fixture Placement



[2] Expander sequence



Review

Kang Yik Je (Director of NY dental surgery)

Although I have used many kits, Haenaem Kit is a convenient kit that can improve bone quality as well as fast speed and stability.

It has already been more than 10 years since the development of various devices for the height adjustment in Korea, and it is thought that it is a product that brings together the advantages of various devices and adds convenience to use.

Jo Jae Beom (Director of Rooted dental surgery)

Many directors complain of discomfort and fear among implant surgery, especially in maxillary sinus lift surgery.

I also used several equipment to solve this. After meeting Haenamber, my surgery became simple and comfortable. It is recommended to try it without worrying.

Kim Si Seok (Director of Rooted dental surgery)

The Osseodensification method makes the Sinus procedure very simple. In particular, I think the biggest advantage is that the Crestal Approach can be operated without burden even for patients with few remaining bones, and that bone densification can be induced without bone loss.

Park Hyoung Mok (Director of Soo San dental surgery)

As a result of performing maxillary sinus lift using the Hanamber kit, bone loss was small and membrane perforation in the maxillary sinus rarely occurred, so membrane lifting was easy without applying external force. It is a product that I would like to recommend because it is possible to perform a lift with only autogenous bone without using different bones.

Jo Seung Heon (Director of Saint dental surgery)

Due to its unique design, Hanamber has the characteristic that bone chips that have been removed during drilling are not removed from the outside, but into the drilling hole. This makes it easier to obtain initial fixation by increasing the bone density by increasing the bone density, or if the maxillary sinus is slightly perforated, the bone chip is inserted into the maxillary sinus during drilling, enabling safe maxillary sinus elevation. In addition, since the drilling is quiet and quiet, it is a great help to maintain the path, and when using other drills, thin bone fragments pop out and the direction of the next drilling or when planting a fixture may change.

On the other hand, It is remarkable for decreases such a risk by using Haenaem Bur. The design of the preparation surface is also important, and the degree of tapering of the drill seems to affect it. Since the first use of Haenaem Bur Kit, the use of implant manufacturers' drills has been significantly reduced.

There is no longer a need to use other maxillary sinus kits.

I recommend you try it out.

Woo Dong Hyup (Director of Boston dental surgery)

The Haenaem Bur Kit maximizes the merits by separating the bone expansion Bur and the maxillary sinus Bur by use, and eliminates mistakes due to rotation direction as a familiar surgeon does not change the implant engine settings through forward drilling. Bone quality is enhanced by bone densification, so even when bone quality is poor, loading time can be accelerated, and autogenous bone transplantation through Crestal Approach enables safe surgery such as less swelling and pain reduction after surgery. In addition, it has a safe bone expansion function through a drill specialized in the narrow bone width of the mandible.

I think this kit is a product that can change the game of existing implant procedures.

Essence Tip

1. During drilling with Haenaem Bur, **“Up & Down” pumping motion is mandatory.**
To maximize the densification effect by creating pressure inside as well as naturally pushing the cut particles inward.
2. Must be perforation of the lower cartilage bone with “L263S”
3. Basic sequence until perforation **(You must proceed in this sequence)**
L18S -> L1824S -> S2228S -> L2632S

We have only 7 kind of stopper sleeve drills as L16S, S18S, L18S, L1628S, L1824S, S2228S, L2632S.

It is indicated with white dot on the kit.

Because the 7 kind drills only using for perforation. After perforation of lower cartilage bone, you can drilling until final size without risk. That is reason why we don't make stopper sleeve to big size drills.

4. Example based on 6mm remaining bone,
L18S (Under 2mm) - L1824S (Under 1mm) – S2228S (Under 1mm) – L2632S (Perforation)
– Depth Gauge
@ If doesn't perforate -> L2632S (+1mm over)
**“It is very important to catch the feeling transmitted to the hand when it is perforate.
Typical feeling: momentarily, the resistance of the force applied to the drill weakens, and the feeling of being pushed easily.”**
5. The synergistic effect is even better when our kit is used in parallel with the existing surgical method used.
EX1)
Perforation with Haenaem Bur -> Sinus lifting with water elevation -> Haenaem Bur as final drill
EX2)
Osteotome with mallet (Under 1mm) -> Haenaem Bur for sinus lifting -> Haenaem Bur as final drill

Essence Tip

6. When using bone powder supplement, the Haenaem Bur Drill (RPM50 without irrigation) allows you to easily push the substances inward.
7. Due to the nature of the drill, the water injected is also naturally sucked into the drilling inside
like cut bone particles, so the watering effect is very good, so the densified bones are not damaged.
8. Excellent the power for keeping path.
No slipping of the drilling towards the weak bone.
Therefore, even in a difficult position (outside), the side wall does not burst and can be safely drilling at the desired location
9. Recommend 600 ~ 1200 RPM.
When perforation (On “L2632S”), we recommend 600 RPM for beginner.
For the other drills, they can use 800 ~ 1200 RPM base on their skill.
10. Must be clockwise rotation.
If they use our drill with counterclockwise, The drill will be high cutting performance.
This is very dangerous, and you need to be careful.
11. It is recommended to use a stopper for first-time users during maxillary Sinus case.
12. It is possible to drill without breaking the side wall on the septum part, such as in cases to be implanted together with the tooth extraction.
13. After 2~5 surgical experiences are accumulated for the first user, most of them combine their technique + Haenaem Bur so that the operation is very convenient and quick. Although little experience is required, The more skilled the operator, the greater the productivity and usability over and over.

#002

HaeNaem Bur Kit for Sinus Lifting



• PILOT DRILL



• SINUS DRILL (600~1200 RPM, 600 RPM when proceed perforation of catilage bone)



- A structure in which 70% of the force is transmitted in the traveling direction and 30% of the force is transmitted in the lateral direction.
- Sinus Autografting / Bone quality improvement possible only by drilling.
- During drilling, the remaining bones and cell lines rise to the maxilla at the same time as the maxillary sinus and lower cartilage are perforated.

Highlight

Easy & Safety maxillary sinus autografting
 Early fixture fixation on general implant placement
 Enhance bone density for poor bone quality through bone condensing
 Pain / Swelling / Recovery Period Reduction

1. Depending upon the implant type and diameter, begin with the narrowest haenaem bur(L1824S) with repeatedly bouncing-pumping motion (RPM600~1200)
2. As the next haenaem bur in the osteotomy, bone will be pushed toward the apical end and will begin to gently lift the membrane and autograft bone.
3. Use the sequential "Zero Bone Loss Drill" with repeatedly bouncing-pumping motion to achieve maximum membrane lift of 3mm and reach final desired width for implant placement.

#003

HaeNaem Bur Kit For **Expander**



1. When drilling with the world-patented bone densification drill "Haenaem bur", No bone loss & overflow occurs.
2. Haenaem bur (Expander Drill) expands and densifies bones at once by drilling.
3. The septum is naturally formed without perforation by drilling.
4. It is very safe and easy to place an implant in a location that requires expansion of the septum by using stopper

Expander Drill



E15S

E1523S

E1826S

E2129S

E2432S

E2735


E3038

E3341

E3644

#004

Bone Expander Kit



1. Increase satisfactory of implant surgery through reducing the pain by minimize the flap open with this system.
2. Done by one drilling for one implant system.
3. Easy to get the path, no bone heat .

SAW Cutting Drill




Using this on the very narrow ridge bone and cutting to the right.



High Cutting Strength and Improved blade can sustain as it is even if it used many times.

Initial Drill (RPM 800-1000)



Use Ø1.8 guide drill to make the first hole for using the bone expander drills.

Square Ratchet Wrench



Using with Bone Expander Drill for make a torque easily.

Bone Expander Drill



Expand the hole by using the bone expander drills.

Square Handle




Using the Square Handle with Bone Expander Drill as adapter for by using hand.

Square Adapter




Using this Adapter for change connecting part to Hand-piece type.

Splitter



This helps to expand ridge bone before using each the bone expander drills.

Direction for Use



1. Use the straight drill (HSD1831) to locate implant to be placed.
2. Use SAW Cutting Drill (HSW70/HSW90) on the very narrow ridge bone and split bone a little bit.
3. To help bone expanding easier, put the splitter(HRS) using malleting inside ridge and hold the handle of the splitter and then move it front and back carefully to expanding.
4. Expand the hole by using the bone expander drill (HBE1826) with the square handle (SH).
5. To help bone expanding easier, put the splitter (HRS) using malleting inside ridge and hold the handle of the splitter and then move it front and back carefully to expanding.
6. Expand the hole by using the bone expander drill (HBE2432/HBE2836/HBE3240/HBE3644) with the square handle (SH) and/or the ratchet wrench (HRW).
7. Repeat ③~⑥ to expand the hole.
8. Place implant.


#005

Double A Guide Kit




1. Each Accurate Guide is equipped with six irrigation holes, making it very easy to irrigate for drilling.
2. There are two guide drills in the kit to make it longer use.
3. Three retention holes can be restored even if the fastening parts of the Accurate Guide and Guide Drill are loosened.

Accurate Guide



1. Can make an accurate guide for location to place an implant.
2. Irrigation hole makes it convenient for drilling to prevent bone heating

Accurate Pin



1. Use for implant placement of two or more at the same time.
2. Easy to get exact positions and path.

Bone Trimmer



1. Easy to organize implant placement position after tooth/ teeth extraction.
2. Easy to clean up the alveolar bone area.

Guide Drill (RPM 500-1200)



1. Use for drilling to implant placement position connecting with the Accurate Guide.
2. Two-step structure.

Lindemann Drill (RPM 500-1200)




Tissue Punch



1. Easy removal of tissues during flapless operation.
2. Precisely deleting only the surrounding tissue with centered on the hole created by the guide drill.

Direction for Use


1. Select the Accurate Guide of the appropriate size by visually checking the interference with the adjacent tooth of the implant placement site and appropriate gap.
2. Fasten up to the first step of Guide Drill in Accurate Guide.
3. Attach the fastened Accurate Guide and Guide Drill to the hand-piece.
4. Hand-piece set to 45 ~ 55 Ncm / 500 ~ 1200RPM
5. Check the position of the alveolar bone to be drilled and gently close the side of the Accurate Guide to the side of the adjacent tooth or Accurate Pin.




6. Use both hands to prevent deviation from the target point and path. (Hold the Accurate Guide with the other hand)



7. Begin drilling with irrigation.




8. Insert Accurate Pin with the same size into the hole created after drilling.



9. Repeat steps 1 through 7 as needed.
10. Be careful not to separate Accurate Guide and Guide Drill during use.

#006

Total Sinus Kit



1. The way of most advanced safe and simple for all of sinus lift operation
2. This is able to provide two ways of the sinus lift operation (Crestal approach / Lateral approach)
3. Excellent in safety, simple operation and visual convenience for the sinus lift


Guide Drill (RPM 800-1000)



GD20

Making a first hole to the point of perforation on cortical bone before main drilling.

Crestal Reamer (RPM 800-1000)



CR28

CR33

CR38

The flattened end tip of the crestal reamer minimizes damage to the membrane.

Lateral Drill (LD65)



LD65

1. Not in case of using window drills, use creating window perforation easier and safety.
2. 2 steps stopper provides more precision perforation to minimize damage to the membrane

Window Drill (RPM 800-1000)



WD65-1

WD65-2

WD65-3

1. When attempting the window perforation of the cartilage of the maxillary sinus, these make it easy to find centering of bone hole which made by the crestal reamer.
2. There are 0.5mm sequential differences (1mm~3mm) that are able to make the window perforation easier.

Bone Carrier (Crestal)



BCC

In case of crestal approach sinus lift, Inserting the bone graft inside of the maxillary sinus with the bone condenser.

Side Cutter



SC2550

Use this instrument, in case of the window cutting surface is not flat or/and rugged.

Aqua Tip



AT3050

To elevate the separated membrane from the maxillary bone, use to connect with silicone tube and syringe for injecting the saline solution.

Bone Carrier (Lateral)



BCL

In case of lateral approach sinus lift, Inserting the bone graft inside of the maxillary sinus with the bone condenser.

Bone Condenser / Depth Gauge



DG0246

1. Measure the elevated depth of the membrane through the band marking and using stoppers.
2. To push the bone graft to inside of the maxillary sinus

Stopper



SSTP020

SSTP030

SSTP040

SSTP050

SSTP060


SSTP070

1. Connecting with a drill to the same length of the cartilage height of maxillary sinus which is measured by CT
2. Connecting with the depth gauge to measure the depth of the elevated membrane


Direction for Use

The Lateral Drill is kind of reamer for perforation with depth guide stopper. The stopper of Lateral drill can be adjust depth level for more safe drilling.

<Lateral Approach Sinus Lift>



<Crestal Approach Sinus Lift>



- 1) Using GD20, designate the location to be drilled on the site where the implant is to be placed.
- 2) Attach a stopper (SSTP 020~070) of an appropriate size according to the thickness of the affected area with the crestal reamer (SRD28, SRD33, SRD38).
- 3) Drill to the position where the membrane of the maxillary sinus is raised by 1mm with the clamped crestal reamer.
- 4) Attach a silicon tube to AT3050 to raise the membrane to an appropriate position.
- 5) Fill the raised space with artificial bones or autogenous bones using BCL and BCC.
- 6) Place the implant.

#007

Crestal Approach Sinus Pro Kit



1. The way of most advanced safe and simple for all of sinus lift operation
2. This is able to provide two ways of the sinus lift operation (Crestal approach / Lateral approach)
3. Excellent in safety, simple operation and visual convenience for the sinus lift

Point Drill (RPM 400-500)



Locate the point of perforation on cortical bone.

PD20

Guide Drill (RPM 800-1200)



Making a first hole to the point of perforation on cortical bone before main drilling.

GD20

Aqua Tip-Elevation



To elevate the separated membrane from the maxillary bone, use to connect with silicone tube and syringe for injecting the saline solution

AT3050

Silicone Tube



SCT

Connecting to the Aqua Tip – Elevation in order to inject the saline solution

Aqua Tip – Dual Action



1. Bone and membrane can be separated easily by omitting the saline solution in both directions at the same time.
2. You can lift the separated membrane directly with the top part made of silicon as well.
3. Designed to make it easy to lift the membrane, it has a world patent.

ATDA

Sinus R Drill (RPM 200-400)



1. It has 6 outer diameters and can be selected according to various clinical cases
2. The rounded tip of drills minimizes damage to the membrane during the operation

Bone Carrier



BCC

Inserting the bone graft inside of the maxillary sinus with the bone pusher

Bone Pusher / Depth Gauge



DG0315


1. Measure the elevated depth of the membrane through the band marking and using stoppers
2. To push the bone graft to inside of the maxillary sinus

Stopper



1. Connecting with a drill to drill to the same length of the cartilage height of maxillary sinus which is measured by CT
2. Connecting with the depth gauge to measure the depth of the elevated membrane

Direction for Use



1. Using PD20 and GD20, designate the drilling position in the area where the implant is to be placed.
2. Connect the Sinus R Drill (SRD28, SRD31, SRD33, SRD36, SRD38, SRD41) to the appropriate size stopper (SSTP020~100) according to the thickness of the affected area.
3. Drill to the position where the membrane of the maxillary sinus is raised by 1mm with the fastened Sinus R Drill.
4. Aqua Tip – Connect the Silicone Tube to Dual Action (ATDA), insert it into the perforated affected area, and inject water to separate the membrane and bone.
5. Aqua Tip – Connect the Silicone Tube to the elevation (AT3050), insert it into the perforated affected area, and inject water to elevate the membrane to an appropriate position.
6. Fill the raised space with artificial bone or autologous bone using BCC and DG0315.
7. Place the implant

#008

One Drilling System Kit



1. Increase satisfactory of implant surgery through reducing the pain by minimize the flap open with this system.
2. Done by one drilling for one implant system.
3. Easy to get the path, no bone heat .
4. Able to collect self-generated bone.

One Drill (RPM 800-1000)



The drill is made ergonomically and provide not only safe drilling but also the size of implant you would like in one drilling instead of drilling many times.

< The benefit of one drill >

1. Able to make the implant size in one drilling.
2. No bone heat.
3. Able to collect self-bone.
4. Big save of the surgery time.

Guide Drill (RPM 800-1000)



GD3248

The first drill you have to use and very good for the setting up the implant position.

Bone Pusher



BP

Can take out the bone inside of the drill.

Bone Remover (Depth Gauge)



BR2116

1. Checking the depth by using line marking on the Bone remover.
2. Removing residue bone on created hole.

Stopper



9 kind of depth you can have.

Opener OUT



OPO

Tools to disjoint the one drill as a shank part and a flute part.

Pincette



PS


Tool to remove remaining bone after using the one drill.

Direction for Use




1. Use the guide drill (GD3248) to be careful on the slippery strong bone surface after minimized flap opening.
2. Use the one drill of the same size drill according to the implant size you would like to place.
3. Check the depth with the bone remover (BR2116). If the depth is not deep enough, use the bone remover to remove the leftover bone with stopper.
4. In case of when you find out very hard bone (D1 or D2 bone), you may use one size bigger drill and put it in half only after using the same size drill of the implant size. In case of when you find out very soft bone, you may use one size smaller drill than the implant size.

How to Remove Bone Chip & Cleaning



#009

Total Remove Kit



Most easy way to remove broken screw & fixture by using Total Remove Kit.

Screw Remove Part

Implant Guide Using screw remove drill/tap/bar depend on situation in reverse side to remove broken screw with guide.



Screw Remove Drill



Screw Remove Tap




Screw Remove Bar




Fixture Remove Part

Fixture Remove




Using fixture remove take fixture and then reverse side to remove fixture by attached ratchet wrench and adapter.



Cover Screw Remove



Using cover screw remover when you find difficulty to open the closing screw.




Instruments

Square Adapter




Hand tool to use fixture remove for handy type.



Square Ratchet Wrench



Hand tool to use fixture remove for high torque.



#010

Solid Screw Kit For GBR


HNGBR

▣ Solid Bone Screw

| Item Code | Diameter | Length (mm) | Color | Q'ty | |
|-----------|----------|-------------|-----------|------|--|
| BC1403 | Ø1.4 | 3 | Green | 5pcs | |
| BC1404 | | 4 | Purple | | |
| BC1406 | | 6 | Yellow | | |
| BC1408 | | 8 | Blue | | |
| BC1603 | Ø1.6 | 3 | Yellow | | |
| BC1604 | | 4 | Dark Blue | | |
| BC1606 | | 6 | Blue | | |
| BC1608 | | 8 | Purple | | |
| BC1610 | | 10 | Green | | |

▣ Solid Bone Tac

| Item Code | Diameter | Length (mm) | Color | Q'ty | |
|-----------|----------|-------------|--------|-------|--|
| BT2535 | Ø2.5 | 3.5 | Blue | 10pcs | |
| BT2545 | | 4.5 | Yellow | | |

▣ Dome Screw


| Item Code | Diameter | Length (mm) | Color | Q'ty | |
|-----------|----------|-------------|--------|------|--|
| DS1511 | Ø5.0 | 11 | Purple | 4pcs | |
| DS1509 | | 9 | Yellow | | |
| DS1507 | | 7 | Blue | | |

▣ Bone Tac Handle



BTH

▣ Bone Tac Holder



BTHD

▣ Pilot Drill




GPD10

▣ Screw Handle



SH

▣ Screw Holder (Hand-piece)



SHH

▣ Screw Holder (Manual)



SHM

#011

Bone Collector Kit



HNBCK

Highlight

1. Choose appropriate bone chip maker drill and stopper in accordance with the size of bone collecting area.
2. Attach the first part of 5mm stopper to the bone chip maker drill.
3. Start drilling with 1cc of irrigation (Recommended 500RPM).
4. Make sure the status of bone collecting with the stopper fully attached.

Bone Collector



BC30

BC37

BC44

BC50

Drill Stopper - 10mm



BCSTP30-10

BCSTP37-10

BCSTP44-10

BCSTP50-10

Drill Stopper - 5mm



BCSTP30-5

BCSTP37-5

BCSTP44-5

BCSTP50-5

#012

V-Bone Collector Set



1. Provides a smooth drilling experience based on excellent cutting effect.
2. It is possible to visually check the amount to be collected with a transparent cap designed with elasticity, and a large amount of bone can be collected easily and conveniently.
3. Easy cleaning and storage with easy attachment and detachment of drill and cap.

Highlight

1. Attach the cap to the drill and fix it to the handpiece.
 2. Drill until the middle part of the cap protrudes (depth about 5mm).
 3. When bone collection is complete, remove the cap and transfer the contents to a separate storage container and repeat the same sequence according to the required amount.
- ※ Cap can be sterilized.

The Components



6.0mm

BC60



VBCC

V-Bone Collector Drill



V-Bone Collector Drill with Cap



After Bone Collecting with Cap



After Bone Collecting without Cap



#013 Initial Kit



HNITK

1. Product that contain only initial drills that can be used for almost implant brands and shapes.
2. Provide 5 stoppers for drilling as an appropriate depth.
3. This can reduce the burden on purchasing the implant surgical kit.

Highlight

1. Point Drill : Using a Point Drill to indicate starting point for placement implant.
2. Initial Drill : This is the drill to expand right after guide drill.
3. Lindemann Drill : This can change the direction of the implant bed and widen the implant bed a little bit.
4. Bone Trimmer: Flattening / Trimming / Removing of hard tissue, tooth and bone.

☛ Bone Trimmer ☛ Point Drill ☛ Lindemann Drill ☛ Initial Drill ☛ Drill Stopper



#014

Bone Mill Kit



HNBMK

1. The guide assembly type drill makes it convenient to combine/disconnect the guide, and it is easy to clean and store.
2. No damage to the conical taper inside the Fixture because the guide support does not rotate during drilling.
3. Solve foundation hole creation and bone mill at once with only drilling using a combination drill

🦷 Bone Mill Guide Narrow



Ø4.0 Ø4.5

HBMG40

HBMG45

If the abutment cannot be properly fixed due to the interference of adjacent bones during the fastening of the abutment after fixture placement, the bone mill guide is drilled at a low speed of less than 100rpm to gradually cut the excess bones around the implantation area.

🦷 Bone Mill Guide Regular



Ø5.0 Ø6.0

HBMG50

HBMG60

If the abutment cannot be properly fixed due to the interference of adjacent bones during the fastening of the abutment after fixture placement, the bone mill guide is drilled at a low speed of less than 100rpm to gradually cut the excess bones around the implantation area.

🦷 Peanut Trimmer



Ø6.0

HPTD60

Peanut Trimmer can be used flexibly without affecting from the location and angle etc.

🦷 Bone Trimmer



Ø5.0

FRB50

Flattening / Trimming / Removing of hard gingiva, tooth and bone.

#015

Tissue Punch Kit



HNTPK

1. A-Type : Before drilling, these can be used to take out the tissue.
2. B-Type : These can be used for same time to take out the tissue and make a 5mm guide hole.
3. C-Type : These can remove the tissue and there is no need extra action to remove residual tissue.

Pusher



PSH

Highlight

1. Choose an appropriate component in accordance with the size of implant placement.
2. When you choose an appropriate component, you can choose B-type if you would like to make a guide hole.
3. Remove tissue with irrigation by using hand-piece

A-Type




HTIPU35

HTIPU45

HTIPU55

B-Type



HTIPU40G

HTIPU50G

C-Type



HTIPU40C

HTIPU50C

#016

Round Bur Kit


HNRBK

The Components




#017

Trephine Kit


HNTHK

Direction of use



You can choose one among its 6 trephine drill in accordance with various usage and diameter.

| ITEM CODE | INNER DIAMETER | ITEM CODE | INNER DIAMETER |
|-----------|----------------|-----------|----------------|
| HTD3015 | Ø2.4 | HTD6015 | Ø5.2 |
| HTD4015 | Ø3.4 | HTD7015 | Ø6.2 |
| HTD5015 | Ø4.2 | HTD8015 | Ø7.2 |







STERI-OS5 IMPLANT Co Ltd

Index

STERI-OS5 IMPLANT Co Ltd, 200, 201, 202, 203, 204, 205, 206, 207, 208, 209, 210, 211, 212, 213, 214, 215, 216, 217, 218, 219, 220, 221, 222, 223, 224, 225, 226, 227, 228, 229, 230, 231, 232, 233, 234, 235, 236, 237, 238, 239, 240, 241, 242, 243, 244, 245, 246, 247, 248, 249, 250, 251, 252, 253, 254, 255, 256, 257, 258, 259, 260, 261, 262, 263, 264, 265, 266, 267, 268, 269, 270, 271, 272, 273, 274, 275, 276, 277, 278, 279, 280, 281, 282, 283, 284, 285, 286, 287, 288, 289, 290, 291, 292, 293, 294, 295, 296, 297, 298, 299, 300, 301, 302, 303, 304, 305, 306, 307, 308, 309, 310, 311, 312, 313, 314, 315, 316, 317, 318, 319, 320, 321, 322, 323, 324, 325, 326, 327, 328, 329, 330, 331, 332, 333, 334, 335, 336, 337, 338, 339, 340, 341, 342, 343, 344, 345, 346, 347, 348, 349, 350, 351, 352, 353, 354, 355, 356, 357, 358, 359, 360, 361, 362, 363, 364, 365, 366, 367, 368, 369, 370, 371, 372, 373, 374, 375, 376, 377, 378, 379, 380, 381, 382, 383, 384, 385, 386, 387, 388, 389, 390, 391, 392, 393, 394, 395, 396, 397, 398, 399, 400, 401, 402, 403, 404, 405, 406, 407, 408, 409, 410, 411, 412, 413, 414, 415, 416, 417, 418, 419, 420, 421, 422, 423, 424, 425, 426, 427, 428, 429, 430, 431, 432, 433, 434, 435, 436, 437, 438, 439, 440, 441, 442, 443, 444, 445, 446, 447, 448, 449, 450, 451, 452, 453, 454, 455, 456, 457, 458, 459, 460, 461, 462, 463, 464, 465, 466, 467, 468, 469, 470, 471, 472, 473, 474, 475, 476, 477, 478, 479, 480, 481, 482, 483, 484, 485, 486, 487, 488, 489, 490, 491, 492, 493, 494, 495, 496, 497, 498, 499, 500, 501, 502, 503, 504, 505, 506, 507, 508, 509, 510, 511, 512, 513, 514, 515, 516, 517, 518, 519, 520, 521, 522, 523, 524, 525, 526, 527, 528, 529, 530, 531, 532, 533, 534, 535, 536, 537, 538, 539, 540, 541, 542, 543, 544, 545, 546, 547, 548, 549, 550, 551, 552, 553, 554, 555, 556, 557, 558, 559, 560, 561, 562, 563, 564, 565, 566, 567, 568, 569, 570, 571, 572, 573, 574, 575, 576, 577, 578, 579, 580, 581, 582, 583, 584, 585, 586, 587, 588, 589, 590, 591, 592, 593, 594, 595, 596, 597, 598, 599, 600, 601, 602, 603, 604, 605, 606, 607, 608, 609, 610, 611, 612, 613, 614, 615, 616, 617, 618, 619, 620, 621, 622, 623, 624, 625, 626, 627, 628, 629, 630, 631, 632, 633, 634, 635, 636, 637, 638, 639, 640, 641, 642, 643, 644, 645, 646, 647, 648, 649, 650, 651, 652, 653, 654, 655, 656, 657, 658, 659, 660, 661, 662, 663, 664, 665, 666, 667, 668, 669, 670, 671, 672, 673, 674, 675, 676, 677, 678, 679, 680, 681, 682, 683, 684, 685, 686, 687, 688, 689, 690, 691, 692, 693, 694, 695, 696, 697, 698, 699, 700, 701, 702, 703, 704, 705, 706, 707, 708, 709, 710, 711, 712, 713, 714, 715, 716, 717, 718, 719, 720, 721, 722, 723, 724, 725, 726, 727, 728, 729, 730, 731, 732, 733, 734, 735, 736, 737, 738, 739, 740, 741, 742, 743, 744, 745, 746, 747, 748, 749, 750, 751, 752, 753, 754, 755, 756, 757, 758, 759, 760, 761, 762, 763, 764, 765, 766, 767, 768, 769, 770, 771, 772, 773, 774, 775, 776, 777, 778, 779, 780, 781, 782, 783, 784, 785, 786, 787, 788, 789, 790, 791, 792, 793, 794, 795, 796, 797, 798, 799, 800, 801, 802, 803, 804, 805, 806, 807, 808, 809, 810, 811, 812, 813, 814, 815, 816, 817, 818, 819, 820, 821, 822, 823, 824, 825, 826, 827, 828, 829, 830, 831, 832, 833, 834, 835, 836, 837, 838, 839, 840, 841, 842, 843, 844, 845, 846, 847, 848, 849, 850, 851, 852, 853, 854, 855, 856, 857, 858, 859, 860, 861, 862, 863, 864, 865, 866, 867, 868, 869, 870, 871, 872, 873, 874, 875, 876, 877, 878, 879, 880, 881, 882, 883, 884, 885, 886, 887, 888, 889, 890, 891, 892, 893, 894, 895, 896, 897, 898, 899, 900, 901, 902, 903, 904, 905, 906, 907, 908, 909, 910, 911, 912, 913, 914, 915, 916, 917, 918, 919, 920, 921, 922, 923, 924, 925, 926, 927, 928, 929, 930, 931, 932, 933, 934, 935, 936, 937, 938, 939, 940, 941, 942, 943, 944, 945, 946, 947, 948, 949, 950, 951, 952, 953, 954, 955, 956, 957, 958, 959, 960, 961, 962, 963, 964, 965, 966, 967, 968, 969, 970, 971, 972, 973, 974, 975, 976, 977, 978, 979, 980, 981, 982, 983, 984, 985, 986, 987, 988, 989, 990, 991, 992, 993, 994, 995, 996, 997, 998, 999, 1000

www.steri-ossimplant.co.uk