


# PROSTHETICS FLOW DIAGRAMS.

Overdenture Restoration



# STERIOSS IMPLANT Reference Code.

## Ball Abutment



- Used in creating stud type overdenture prosthetics
- Gold coloring for improved aesthetics
- Rotator Abutment driver
- Recommended tightening torque : 30Ncm

| Cuff | Ø | 4.0      |
|------|---|----------|
| 0    |   | AM2BAR00 |
| 1    |   | AM2BAR01 |
| 2    |   | AM2BAR02 |
| 3    |   | AM2BAR03 |
| 4    |   | AM2BAR04 |
| 5    |   | AM2BAR05 |
| 6    |   | AM2BAR06 |

## Ball Cap



- Used in creating stud type overdenture prosthetics

| Ø | 4.0   |
|---|-------|
|   | AM2BC |

## Ball Lab Analog



| Ø | 3.5   |
|---|-------|
|   | BLA35 |