



PROSTHETICS FLOW DIAGRAMS.

Overdenture Restoration



STERIOSS IMPLANT Reference Code.

Rotator Abutment



- Used in creating Rotator type overdenture prosthetics
- Gold coloring for improved aesthetics
- Rotator Abutment driver
- Recommended tightening torque : 30Ncm

Ø	
Cuff	
1	LOGSR3710
2	LOGSR3720
3	LOGSR3730
4	LOGSR3740
5	LOGSR3750
6	LOGSR3760
7	LOGSR3770
8	LOGSR3780

Impression Cap



- Used in creating Rotator type overdenture prosthetics

Ø	
4.0	ROTIP

Rotator Analog



Ø	
4.0	ROTAL 40

- Achieves Rotator Abutment of the oral cavity on a working model