



PROSTHETICS FLOW DIAGRAMS.

Overdenture & Cement Retained Restoration



STERIOSS IMPLANT Reference Code.

Rotator Abutment



Ø	
Cuff	
0	LGSSSR4800
1	LGSSSR4810
2	LGSSSR4820
3	LGSSSR4830
4	LGSSSR4840

- Used in creating stud type overdenture prosthetics
- Gold coloring for improved aesthetics
- Ball Abutment driver
- Recommended tightening torque : 30Ncm

Impression Cap



Ø	4.0
	ROTIP

- Used in creating Rotator type overdenture prosthetics

Rotator Lab Analog



Ø	4.0
	ROTAL 40

- Achieves Rotator Abutment of the oral cavity on a working model